



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Rasotherm beaker low form with spout, (Boro 3.3), 10 ml

Special Price

**€2.10** was

**€2.72** As low as

**€1.99**

**plus VAT &  
Shipping**

### Product Images



## Description

---

The RASOTHERM® **beakers**, low form with spout are made of especially heat & chemical resistant borosilicate glass 3.3. With graduation.

All beakers meet the highest quality standards and are manufactured according to the specifications of the current ISO standard 3819.

Thus, the beakers are suitable for many purposes in your laboratory.

## Additional Information

---

No.	1-0150
Manufacturer (Brand)	Rasotherm
EAN	4058072080846
autoclavable at 121°C	Yes
diameter	26 mm
Color	transparent (clear)
sterile	No
Form	low form
Height	35 mm
Cups and scales properties	with spout suitable for food
Material	Borosilicate glass
Material cup	Borosilicate glass
Type cup	Beakers
Volume cup	10 ml

