



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Rasotherm beaker low form with spout, (Boro 3.3), 5000 ml

Special Price

**€29.40** was

**€37.70**

**plus VAT &  
Shipping**

## Product Images

---



## Description

---

The RASOTHERM® **beakers**, low form with spout are made of especially heat & chemical resistant borosilicate glass 3.3. With graduation.

All beakers meet the highest quality standards and are manufactured according to the specifications of the current ISO standard 3819.

Thus, the beakers are suitable for many purposes in your laboratory.

## Additional Information

---

No.	1-0162
Manufacturer (Brand)	Rasotherm
EAN	4058072080969
autoclavable at 121°C	Yes
diameter	170 mm
Color	transparent (clear)
sterile	No
Form	low form
Height	270 mm
Cups and scales properties	with spout graduated suitable for foodstuffs
Material	Borosilicate glass
Color graduation	white
Material cup	Borosilicate glass
Type cup	Beakers
Volume cup	5000 ml

