



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Rasotherm beaker high form with spout, (Boro 3.3), 800 ml

Special Price

€5.40 was

€6.93 As low as

€4.94

**plus VAT &
Shipping**



Description

The RASOTHERM® **beakers**, high form with spout are made of especially heat & chemical resistant borosilicate glass 3.3.

With graduation.

All beakers meet the highest quality standards and are manufactured according to the specifications of the current ISO standard 3819.

Thus, the beakers are suitable for many purposes in your laboratory.

Additional Information

No.	1-0170
Manufacturer (Brand)	Rasotherm
EAN	4058072081041
autoclavable at 121°C	Yes
diameter	90 mm
Color	transparent (clear)
sterile	No
Form	high shape
Height	175 mm
Cups and scales properties	with spout graduated suitable for foodstuffs
Material	Borosilicate glass
Color graduation	white
Material cup	Borosilicate glass
Type cup	Beakers
Volume cup	800 ml

