



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Disposable paper cups 500 ml, PE-coated, 100 pcs./pack

€57.00
plus VAT &
Shipping

Product Images



Description

Indispensable for filling, measuring, assessing or mixing liquids or powders.

The PE coating prevents the cup from becoming wet or damp.

Polyethylene behaves groundwater-neutral and is compostable. Burning produces only water vapor and harmless CO₂, food grade.

Dimensions:

Order No. 1-7103

157 mm high, 90 mm Ø

Additional Information

No.	1-7103
Manufacturer (Brand)	neoLab
EAN	4058072045043
autoclavable at 121°C	No
diameter	90 mm
Color	White
sterile	No
Height	157 mm
Cups and scales properties	for 1xuse
Material	Paper PE coated
Material cup	Paper PE coated
Type cup	1xMug
Volume cup	500 ml

