



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® beaker 100 ml, screw cap yellow,
single packed, 150 pcs./pack**

**€90.50
plus VAT &
Shipping**

Product Images



Description

Transparent, graduated beakers made of PP with 100 ml capacity.

With yellow screw cap, leak-proof and liquid-tight, made of LD-PE (sterile).

Pressure-resistant and unbreakable, stable and sturdy.

Large opening for easy and safe filling.

Writing area on cup and frosted labeling area on lid.

Additional Information

No.	2-0102
Manufacturer (Brand)	neoLab
EAN	4058072047429
autoclavable at 121°C	No
diameter	56 mm
Color	transparent (milky)
Material	Polypropylene (PP)
sterile	Yes
Height	70 mm
Cups and scales properties	with raised graduation with closure/lid with labeling field for lxuse
Type closure/ lid	Screw cap
TBST MAX	100 °C
Material	Polypropylene (PP)
Type of sterility	radiation sterile
Color closure/ lid	yellow
Format/ Packing	individually packed in carton/ dispenser
Material cup	Polypropylene (PP)
Material closure/ lid	Low-density polyethylene (LDPE)
Type cup	Mug with lid
Volume cup	100 ml

