



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Beakers 100 ml with screw cap yellow,  
assembled, non-sterile, 48x5 (=240) pcs./pack**

**€96.10  
plus VAT &  
Shipping**

## Product Images



## Description

---

Transparent, graduated cups made of PP with 100 ml capacity.

With yellow screw cap, leak-proof and liquid-tight, made of LD-PE.

Pressure-resistant and unbreakable, stable and sturdy.

Large opening for easy and safe filling.

Writing area on cup and frosted labeling area on lid.

## Additional Information

---

No.	2-0101
Manufacturer (Brand)	neoLab
EAN	4058072047412
autoclavable at 121°C	No
diameter	56 mm
Color	transparent (milky)
Material	Polypropylene (PP)
sterile	No
Height	70 mm
Cups and scales properties	with raised graduation with closure/lid with labeling field for 1xuse
Type closure/ lid	Screw cap
TBST MAX	100 °C
Material	Polypropylene (PP)
Color closure/ lid	yellow
Format/ Packing	individually packed in carton/ dispenser
Material cup	Polypropylene (PP)
Material closure/ lid	Low-density polyethylene (LDPE)
Type cup	Mug with lid
Volume cup	100 ml

