



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Multi-purpose beaker with screw cap, 60 ml, PP, non-sterile

**€77.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Multipurpose beaker with screw cap made of polypropylene (PP) or polystyrene (PS). Suitable for sampling, transport and storage of solids, pastes and liquids. Leak-proof with red screw cap non-sterile or individually sterile packed

## Additional Information

---

No.	2-0162
Manufacturer (Brand)	neoLab
EAN	4058072047702
autoclavable at 121°C	No
diameter	42,5 mm
Color	transparent (milky)
Material	Polypropylene (PP)
sterile	No
Height	65.5 mm
Cups and scales properties	with closure/lid with labeling field graduated for 1xuse spill-proof
Type closure/ lid	Screw cap
TBST MAX	100 °C
Material	Polypropylene (PP)
Figure	2
Color graduation	black
Color closure/ lid	red
Format/ Packing	single. in carton/ dispenser
Material cup	Polypropylene (PP)
Material closure/ lid	Polypropylene (PP)
Type cup	Mug with lid
Volume cup	60 ml

