



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Multi-purpose beaker with screw cap, 120 ml, PP, sterile

€80.10
plus VAT &
Shipping

Product Images



Description

Multipurpose beaker with screw cap made of polypropylene (PP) or polystyrene (PS). Suitable for sampling, transport and storage of solids, pastes and liquids. Leak-proof with red screw cap non-sterile or individually sterile packed

Additional Information

| | |
|----------------------------|--|
| No. | 2-0165 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072047733 |
| autoclavable at 121°C | No |
| diameter | 62,5 mm |
| Color | transparent (milky) |
| Material | Polypropylene (PP) |
| sterile | Yes |
| Height | 72.5 mm |
| Cups and scales properties | with closure/lid with labeling field graduated for 1xuse spill-proof |
| Type closure/ lid | Screw cap |
| TBST MAX | 100 °C |
| Material | Polypropylene (PP) |
| Figure | 3 |
| Color graduation | black |
| Type of sterility | machine sterile |
| Color closure/ lid | red |
| Format/ Packing | individually packed |
| Material cup | Polypropylene (PP) |
| Material closure/ lid | Polypropylene (PP) |
| Type cup | Mug with lid |
| Volume cup | 120 ml |

