



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Multi-purpose beaker with screw cap, 120 ml, PP, non-sterile

€65.40
plus VAT &
Shipping

Product Images



Description

Multipurpose beaker with screw cap made of polypropylene (PP) or polystyrene (PS). Suitable for sampling, transport and storage of solids, pastes and liquids. Leak-proof with red screw cap non-sterile or individually sterile packed

Additional Information

No.	2-0164
Manufacturer (Brand)	neoLab
EAN	4058072047726
autoclavable at 121°C	No
diameter	62,5 mm
Color	transparent (milky)
Material	Polypropylene (PP)
sterile	No
Height	72.5 mm
Cups and scales properties	with closure/lid with labeling field graduated for 1xuse spill-proof
Type closure/ lid	Screw cap
TBST MAX	100 °C
Material	Polypropylene (PP)
Figure	3
Color graduation	black
Color closure/ lid	red
Format/ Packing	single. in carton/ dispenser
Material cup	Polypropylene (PP)
Material closure/ lid	Polypropylene (PP)
Type cup	Mug with lid
Volume cup	120 ml

