



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® disposable beakers made of PS, 5 ml,  
1000 pcs./pack**

**€71.50  
plus VAT &  
Shipping**

## Product Images



## Description

---

Sample cup made of PS.  
Operating range from -20°C to +70°C.

## Additional Information

---

No.	2-1555
Manufacturer (Brand)	neoLab
EAN	4058072057381
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	32 mm
Color	transparent (milky)
sterile	No
Height	18 mm
Cups and scales properties	with pouring edge
TBST MAX	70 °C
TBST MIN	-10 °C
Material	Polystyrene (PS)
Material cup	Polystyrene (PS)
Type cup	1xMug
Volume cup	5 ml

