



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Disposable beakers made of PS, 150 ml, 100 pcs./pack

**€76.80
plus VAT &
Shipping**

Product Images



Description

Sample cup made of PS.
Operating range from -20°C to +70°C.

Additional Information

No.	2-1559
Manufacturer (Brand)	neoLab
EAN	4058072057428
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	82 mm
Color	transparent (milky)
sterile	No
Height	65 mm
Cups and scales properties	with pouring edge
TBST MAX	70 °C
TBST MIN	-10 °C
Material	Polystyrene (PS)
Material cup	Polystyrene (PS)
Type cup	1xMug
Volume cup	150 ml

