



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® beakers made of PP, 30 ml, yellow, 75
pcs./pack**

**€9.17
plus VAT &
Shipping**

Product Images



Description

Cup made of PP, highly transparent, the content is optimally visible.

Pressure-resistant and unbreakable.

Content: 30 ml

Accessories: Snap-on lid made of PP (Order No. 2-0114)

Additional Information

No.	2-0113
Manufacturer (Brand)	neoLab
EAN	4058072047474
autoclavable at 121°C	Yes
diameter	35 mm
Color	yellow
sterile	No
Height	40 mm
Cups and scales properties	with special chemical resistance
TBST MAX	120 °C
TBST MIN	-10 °C
Material	Polypropylene (PP)
Material cup	Polypropylene (PP)
Type cup	Mug with lid
Volume cup	30 ml

