



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Sample storage jar, with rubber membrane, 60 ml

Special Price

€4.50 was

€6.75

**plus VAT &
Shipping**

Product Images



Description

The sample storage jars have an integrated rubber membrane in the lid that can be easily pierced with a needle. The jars are easy to open and close with one hand. They are made of high purity polypropylene and are excellent for use with food. They were manufactured under clean room conditions. They comply with CE standard 1935/2004.

Additional Information

No.	2-5264
Manufacturer (Brand)	neoLab
EAN	4058072182755
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	40 mm
Color	transparent (milky)
Material	Acrylonitrile butadiene styrene (ABS)
sterile	No
Height	75 mm
Cups and scales properties	with closure/lid with labeling field suitable for food for 1xuse
Type closure/ lid	Slip closure
Material	Acrylonitrile butadiene styrene (ABS)
Color closure/ lid	transparent (crystal clear)
Format/ Packing	in bag
Material cup	Acrylonitrile butadiene styrene (ABS)
Material closure/ lid	Acrylonitrile butadiene styrene (ABS)
Type cup	Mug with lid
Volume cup	60 ml

