



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Moonlab® plastic beaker w. spout o. handle, 2000 ml, graduated, natural, PP, 1 pc

€6.50
plus VAT &
Shipping

Product Images



Description

Plastic PP cup with spout, without handle

This transparent sturdy cup made of polypropylene is the cheap alternative to glass. It is chemically resistant to many acids, alkalis and solvents.

- with spout
- graduated
- autoclavable

Additional Information

| | |
|----------------------------|--|
| No. | 4-0043 |
| Manufacturer (Brand) | Moonlab |
| EAN | 4058072085483 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| diameter | 130 mm |
| Color | transparent (milky) |
| sterile | No |
| Form | around |
| Height | 192 mm |
| Cups and scales properties | with spout graduated suitable for foodstuffs |
| Material | Polypropylene (PP) |
| Color graduation | black |
| Division | 20 ml |
| Material cup | Polypropylene (PP) |
| Type cup | Measuring cup |
| Volume cup | 2000 ml |

