



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Moonlab® Plastic beaker w. spout o. handle, 100 ml, graduated, natural, PP, 1 pc

**€1.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Plastic PP cup with spout, without handle

This transparent sturdy cup made of polypropylene is the cheap alternative to glass. It is chemically resistant to many acids, alkalis and solvents.

- with spout
- graduated
- autoclavable

## Additional Information

---

No.	4-0039
Manufacturer (Brand)	Moonlab
EAN	4058072085490
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	55 mm
Color	transparent (milky)
sterile	No
Form	around
Height	73 mm
Cups and scales properties	with spout graduated suitable for foodstuffs
Material	Polypropylene (PP)
Color graduation	black
Division	2 ml
Material cup	Polypropylene (PP)
Type cup	Measuring cup
Volume cup	100 ml

