



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Moonlab® plastic beaker w. spout & handle, 5000 ml, graduated, natural, PP, 1 pc

**€6.66**  
**plus VAT & Shipping**

### Product Images



## Description

---

### Plastic cup made of PP with spout and handle

The cheap alternative to glass. They are chemically resistant to many acids and alkalis and some solvents.

- with spout
- graduated
- autoclavable
- with practical handle

## Additional Information

---

No.	4-0051
Manufacturer (Brand)	Moonlab
EAN	4058072085421
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	175 mm
Color	transparent (milky)
sterile	No
Form	around
Height	263 mm
Cups and scales properties	with handle with spout graduated food-safe
Material	Polypropylene (PP)
Color graduation	black
Division	50 ml
Material cup	Polypropylene (PP)
Type cup	Measuring cup
Volume cup	5000 ml

