



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Moonlab® plastic beaker w. spout o. handle, 50 ml, graduated, natural, PP, 1 pc

€1.40
plus VAT &
Shipping

Product Images



Description

Plastic PP cup with spout, without handle

This transparent sturdy cup made of polypropylene is the cheap alternative to glass. It is chemically resistant to many acids, alkalis and solvents.

- with spout
- graduated
- autoclavable

Additional Information

No.	4-0038
Manufacturer (Brand)	Moonlab
EAN	4058072085438
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	47 mm
Color	transparent (milky)
sterile	No
Form	around
Height	59 mm
Cups and scales properties	with spout graduated suitable for foodstuffs
Material	Polypropylene (PP)
Color graduation	black
Material cup	Polypropylene (PP)
Division	2 ml
Type cup	Measuring cup
Volume cup	50 ml

