



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Tripour beaker made of PP, 1000 ml,
pack a 100 pcs.**

**€161.30
plus VAT &
Shipping**

Product Images



Description

Shatterproof cups with three specially shaped spout edges to prevent dripping- thus more safety when pouring the contents.

- translucent
- autoclavable
- largely acid-resistant
- reusable

Material: PP

Additional Information

No.	7-4729
Manufacturer (Brand)	neoLab
EAN	4058072041137
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	115 mm
Color	transparent (clear)
sterile	No
Height	145 mm
Cups and scales properties	graduated
Material	Polypropylene (PP)
Color graduation	transparent (milky)
Division	50 ml
Material cup	Polypropylene (PP)
Type cup	Tripour mug
Volume cup	1000 ml

