




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Tripour beaker made of PP, 400 ml,
pack a 100 pcs.**

**€102.70
plus VAT &
Shipping**

Product Images



Description

Shatterproof cups with three specially shaped spout edges to prevent dripping- thus more safety when pouring the contents.

- translucent
- autoclavable
- largely acid-resistant
- reusable

Material: PP

Additional Information

No.	7-4727
Manufacturer (Brand)	neolab
EAN	4058072041113
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	76 mm
Color	transparent (clear)
sterile	No
Height	108 mm
Cups and scales properties	graduated
Material	Polypropylene (PP)
Color graduation	transparent (milky)
Division	20 ml
Material cup	Polypropylene (PP)
Type cup	Tripour mug
Volume cup	400 ml

