



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Tripour beaker made of PP, 250 ml,  
pack a 100 pcs.**

**€113.60  
plus VAT &  
Shipping**

## Product Images



## Description

---

Shatterproof cups with three specially shaped spout edges to prevent dripping- thus more safety when pouring the contents.

- translucent
- autoclavable
- largely acid-resistant
- reusable

**Material:** PP

## Additional Information

---

No.	7-4726
Manufacturer (Brand)	neoLab
EAN	4058072041106
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	76 mm
Color	transparent (clear)
sterile	No
Height	90 mm
Cups and scales properties	graduated
Material	Polypropylene (PP)
Color graduation	transparent (milky)
Division	10 ml
Material cup	Polypropylene (PP)
Type cup	Tripour mug
Volume cup	250 ml

