



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® drawstring bag (PE), transparent, 15 x
22 cm, 100 pcs./pack**

€24.90
**plus VAT &
Shipping**

Product Images



Description

Transparent bags made of polyethylene with white text field and practical cord closure.
Film thickness: 0.05 mm

Additional Information

No.	1-7082
Manufacturer (Brand)	neoLab
EAN	4058072044848
autoclavable at 121°C	No
Color	transparent (clear)
Material	Low-density polyethylene (LDPE)
sterile	No
Bag features	with closure/lid
Wide	150 mm
Length	220 mm
Type bag	Storage bag sample bag
Type closure/ lid	Cord closure
Thickness	0.05 mm

