



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Plastic bags made of LDPE, black, 70 x
100 mm, 100 pcs./pack**

**€18.20
plus VAT &
Shipping**

Product Images



Description

Black LDPE bags for light-sensitive samples, with practical quick closure and labeling field.
By squeezing, the bag is closed lightning fast and tight.
Film thickness: 0.05 mm

Additional Information

No.	1-7130
Manufacturer (Brand)	neoLab
EAN	4058072045241
autoclavable at 121°C	No
Color	Black
Material	Low-density polyethylene (LDPE)
sterile	No
Bag features	with closure/lid with labeling field for light-sensitive specimens
Wide	70 mm
Length	100 mm
Type bag	Sample Bag Storage Bag
Type closure/ lid	Self-adhesive closure
Thickness	0.05 mm

