



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Plastic bags made of LDPE, black, 160 x  
220 mm, 100 pcs./pack**

**€31.50  
plus VAT &  
Shipping**

## Product Images



## Description

---

Black LDPE bags for light-sensitive samples, with practical quick closure and labeling field.  
By squeezing, the bag is closed lightning fast and tight.  
Film thickness: 0.05 mm

## Additional Information

---

No.	1-7133
Manufacturer (Brand)	neoLab
EAN	4058072045272
autoclavable at 121°C	No
Color	Black
Material	Low-density polyethylene (LDPE)
sterile	No
Bag features	with closure/lid with labeling field for light-sensitive specimens
Wide	160 mm
Length	200 mm
Type bag	Sample Bag Storage Bag
Type closure/ lid	Self-adhesive closure
Thickness	0.05 mm

