



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Plastic bag with labeling field, 70 x 100 mm, 100 pcs./pack

**€13.60
plus VAT &
Shipping**

Product Images



Description

LDPE bags with practical quick closure.
Closed in a flash and tight by squeezing the bags together.
Suitable for storing or shipping samples.
Film thickness: 0.05 mm

Additional Information

No.	1-7090
Manufacturer (Brand)	neoLab
EAN	4058072044923
autoclavable at 121°C	No
Color	transparent (clear)
Material	Low-density polyethylene (LDPE)
sterile	No
Bag features	with closure/lid with labeling field
Wide	70 mm
Length	100 mm
Type bag	Storage bag sample bag shipping bag
Type closure/ lid	Pressure cap
Thickness	0.05 mm

