



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Plastic bags 180 x 250 mm, 100
pcs./pack**

**€20.80
plus VAT &
Shipping**

Product Images



Description

LDPE bags with practical quick closure.
Closed in a flash and tight by squeezing the bags together.
Suitable for storing or shipping samples.
Film thickness: 0.05 mm

Additional Information

No.	1-7124
Manufacturer (Brand)	neoLab
EAN	4058072045210
autoclavable at 121°C	No
Color	transparent (clear)
Material	Low-density polyethylene (LDPE)
sterile	No
Bag features	with closure/lid
Wide	180 mm
Length	250 mm
Type bag	Storage bag sample bag shipping bag
Type closure/ lid	Pressure cap
Thickness	0.05 mm

