



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® sealed edge pouches made of
aluminum compound foil, 100 x 250 mm, 100
pcs./pack**

**€107.80
plus VAT &
Shipping**

Product Images



Description

Order No. 7-0080

Sealed edge pouch made of aluminum composite film, 117µm thick, 100 x 250 mm.
Other sizes on request.

We will be pleased to offer you the suitable film sealing device on request.

Additional Information

No.	7-0080
Manufacturer (Brand)	neoLab
EAN	4058072026929
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Aluminum
sterile	No
Bag features	with closure/lid for light-sensitive specimens
Wide	100 mm
Length	250 mm
Type bag	Storage bag
Type closure/ lid	Seal edge closure
Thickness	0.12 mm

