



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® sealed edge pouches made of  
aluminum compound foil, 150 x 300 mm, 100  
pcs./pack**

**€108.80  
plus VAT &  
Shipping**

## Product Images



## Description

---

Order No. 7-0081

Sealed edge bag made of aluminum composite film, 117µm thick, 150 x 300 mm.  
Other sizes on request.

We will be pleased to offer you the suitable film sealing device on request.

## Additional Information

---

No.	7-0081
Manufacturer (Brand)	neoLab
EAN	4058072026936
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Aluminum
sterile	No
Bag features	with closure/lid for light-sensitive specimens
Wide	150 mm
Length	300 mm
Type bag	Storage bag
Type closure/ lid	Seal edge closure
Thickness	0.12 mm

