



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Ziplock bag, closure yellow, size
400x500 mm, 2 x 30 pcs./pack**

**€25.90
plus VAT &
Shipping**

Product Images



Description

Transparent plastic bag with yellow zip closure (sliding closure). Samples can be easily differentiated by the colored closure. Quick and easy to close. Ideal bag for small parts and samples.

- Material LDPE
- Use temperature: -50 to 75°C
- Film thickness from 0.05 to 0.10 mm
- Multiple use

Additional Information

No.	1-7341
Manufacturer (Brand)	neoLab
EAN	4058072046415
autoclavable at 121°C	No
Color	transparent (clear) yellow
Material	Low-density polyethylene (LDPE)
sterile	No
Bag features	with closure/lid
Wide	500 mm
Length	400 mm
Type bag	Storage bag Shipping bag
Type closure/ lid	Duck closure with slider (e.g. Ziplock®)
TBST MAX	75 °C
TBST MIN	-50 °C
Thickness	0.1 mm
Color closure/ lid	yellow

