



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Bag for BagMixer 2-1080, 10 x 15 cm, 50 pcs./pack

Special Price

**€3.10** was

**€4.65**

**plus VAT &  
Shipping**

### Product Images



## Description

---

For extracting, mixing, homogenizing and dispersing. Samples are enclosed in a sterile plastic bag (possibly with the addition of diluent) and clamped in the instrument. The beater blades grind the sample without the sample coming into contact with the mixer. A duration of 30-60 seconds is usually sufficient. The beating speed is adjustable (6-10 beats/second). Housing and mixing chamber made of stainless steel easy to clean due to fully opening closing flap double safety flap with viewing window, watertight adjustable mixing blades for adaptation to sample size and type (patented) mixing blades exchangeable with drip tray Technical data model MiniMix® Timer: 10-360 seconds Dimensions: 35 x 20 x 17.5 cm Weight: 8 kg Matching sterile bags 2-1081 in size 10 x 15 cm.

## Additional Information

---

No.	2-1081
Manufacturer (Brand)	neolab
EAN	4058072055424
Color	transparent (clear)
Material	Low-density polyethylene (LDPE. Multilayer®)
sterile	Yes
Wide	15 mm
Length	10 mm
Type bag	Sterile bags for Bagmixer® MiniMix®
for part no.	2-1080

