



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Whirl Pak® plastic bags with title block, self-standing, sterile, PE**

**€71.30**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

These **self-standing Whirl Pak® bags** make sample handling easier. They no longer need to be held, allowing the user to perform important laboratory functions with both hands. The bottom of the self-standing pouches unfolds when a sample is placed inside the pouch. This allows the bag to stand upright independently. Whirl Pak self-standing pouches function just like a bottle or beaker. However, they cost less, are sterilized, shatterproof, and use only a fraction of the storage space that bottles and beakers would require.

## Additional Information

---

No.	1-7284
Manufacturer (Brand)	Whirl Pak
EAN	4058072155131
autoclavable at 121°C	No
Color	transparent (clear)
Material	Low-density polyethylene (LDPE)
Volume	1242 ml
sterile	Yes
Bag features	with closure/lid with labeling field
Wide	150 mm
Length	380 mm
Type bag	Sample Bag Storage Bag
Type closure/ lid	Swing stopper
TBST MAX	82 °C
TBST MIN	-210 °C
Thickness	1.02 mm

