



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Whirl Pak® plastic bags with title block and filter, 300x190 mm, sterile, PE, 50 pcs./pack

**€163.40
plus VAT &
Shipping**

Product Images



Description

Whirl-Pak® bags made of PE with labeling field and filter

The bags have a filter layer made of finely perforated polyethylene to separate the liquid and solids so that the sample can be pipetted easily. The filter is enclosed within the perimeter of the bag on the sides and bottom so that the sample can be placed in the bag on either side. The liquid is transferred through the filter, but the solids remain on one side.

Additional Information

No.	1-7153
Manufacturer (Brand)	Whirl Pak
EAN	4058072045463
autoclavable at 121°C	No
Color	transparent (clear)
Material	Low-density polyethylene (LDPE)
Volume	1650 ml
sterile	Yes
Bag features	with closure/lid with filter with labeling field
Wide	190 mm
Length	300 mm
Type bag	Sample Bag Storage Bag
Type closure/ lid	Swing stopper
TBST MAX	82 °C
TBST MIN	-210 °C
Thickness	1.02 mm

