



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Whirl Pak® plastic bag with title block and filter, 254x380 mm, sterile


Special Price

**€99.00** was

**€148.50**

plus VAT &  
Shipping

### Product Images



## Description

---

### **Whirl-Pak® bags made of PE with labeling field and filter**

The bags have a filter layer made of finely perforated polyethylene to separate the liquid and solids so that the sample can be pipetted easily. The filter is enclosed within the perimeter of the bag on the sides and bottom so that the sample can be placed in the bag on either side. The liquid is transferred through the filter, but the solids remain on one side.

## Additional Information

---

No.	1-7160
Manufacturer (Brand)	Whirl Pak
EAN	4058072156275
autoclavable at 121°C	No
Color	transparent (clear)
Material	Low-density polyethylene (LDPE)
Volume	2721 ml
sterile	Yes
Bag features	with closure/lid with filter with labeling field
Wide	254 mm
Length	380 mm
Type bag	Sample Bag Storage Bag
Type closure/ lid	Swing stopper
TBST MAX	82 °C
TBST MIN	-210 °C
Thickness	1.02 mm

