



neolab Migge GmbH

Rischerstr. 7-9

69123 Heidelberg

Deutschland

+49 (0)6221 / 8442-44

<https://www.neolab.de>

Umsatzsteuer-Identifikationsnummer:


DE 143 450 657



## **Whirl Pak® plastic bag with title block and filter, 375x375 mm, sterile**

€367.20

plus VAT & Shipping



## Description

### Whirl-Pak® bags made of PE with labeling field and filter

The bags have a filter layer made of finely perforated polyethylene to separate the liquid and solids so that the sample can be pipetted easily. The filter is enclosed within the perimeter of the bag on the sides and bottom so that the sample can be placed in the bag on either side. The liquid is transferred through the filter, but the solids remain on one side.

## Additional Information

No.	1-7161
Manufacturer (Brand)	Whirl Pak
EAN	4058072155001
autoclavable at 121°C	No
Color	transparent (clear)
Material	Low-density polyethylene (LDPE)
Volume	4080 ml
sterile	Yes
Bag features	with closure/lid with filter with labeling field
Wide	375 mm
Length	375 mm
Type bag	Sample Bag Storage Bag
Type closure/ lid	Swing stopper
TBST MAX	82 °C
TBST MIN	-210 °C
Thickness	1.02 mm

