



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Whirl Pak® plastic bag with title block, 23 x 11.5
cm (L x W)**

**€40.30
plus VAT &
Shipping**

Product Images



Description

Versatile **sample bags** especially suitable for liquid samples for food, water and wastewater analysis, but also for biological or chemical samples.

To open and close the bag, the filling opening is backed on both sides with a metal clip, which also keeps the opening taut at all times.

For hermetic sealing, hold the bag by the metal clip and spin it around its own axis several times, then fold in the two ends of the metal clip and press tightly. The Whirl-Pak® can be opened and closed as often as required.

A strong weld seam on all sides guarantees additional safety.

Additional Information

No.	1-7140
Manufacturer (Brand)	Whirl Pak
EAN	4058072045333
autoclavable at 121°C	No
Color	transparent (clear)
Material	Low-density polyethylene (LDPE)
Volume	540 ml
sterile	Yes
Bag features	with closure/lid with labeling field
Wide	115 mm
Length	230 mm
Type bag	Sample Bag Storage Bag
Type closure/ lid	Swing stopper
TBST MAX	82 °C
TBST MIN	-210 °C
Thickness	0.06 mm

