



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Whirl Pak® plastic bag, 38 x 15 cm (L x W)

€67.50
plus VAT &
Shipping

Product Images



Description

Versatile **sample bags** especially suitable for liquid samples for food, water and wastewater analysis, but also for biological or chemical samples.

To open and close the bag, the filling opening is backed on both sides with a metal clip which also keeps the opening taut at all times.

For hermetic sealing, hold the bag by the metal clip and spin it around its own axis several times, then fold in the two ends of the metal clip and press tightly. The Whirl-Pak® can be opened and closed as often as required.

A strong weld seam on all sides guarantees additional safety.

Additional Information

No.	1-7142
Manufacturer (Brand)	Whirl Pak
EAN	4058072045357
autoclavable at 121°C	No
Color	transparent (clear)
Material	Low-density polyethylene (LDPE)
Volume	1260 ml
sterile	Yes
Bag features	with closure/lid
Wide	150 mm
Length	380 mm
Type bag	Sample Bag Storage Bag
Type closure/ lid	Swing stopper
TBST MAX	82 °C
TBST MIN	-210 °C
Thickness	0.74 mm

