



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Whirl Pak® plastic bag, 18.5 x 7.5 cm (L x W)

**€29.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Versatile **sample bags** especially suitable for liquid samples for food, water and wastewater analysis, but also for biological or chemical samples.

To open and close the bag, the filling opening is backed on both sides with a metal clip which also keeps the opening taut at all times.

For hermetic sealing, hold the bag by the metal clip and spin it around its own axis several times, then fold in the two ends of the metal clip and press tightly. The Whirl-Pak® can be opened and closed as often as required.

A strong weld seam on all sides guarantees additional safety.

## Additional Information

---

No.	1-7144
Manufacturer (Brand)	Whirl Pak
EAN	4058072045371
autoclavable at 121°C	No
Color	transparent (clear)
Material	Low-density polyethylene (LDPE)
Volume	120 ml
sterile	Yes
Bag features	with closure/lid
Wide	75 mm
Length	185 mm
Type bag	Sample Bag Storage Bag
Type closure/ lid	Swing stopper
TBST MAX	82 °C
TBST MIN	-210 °C
Thickness	0.57 mm

