



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Whirl Pak® plastic pouch with title block, 150 x 230 mm (W x L), sterile

€70.60
plus VAT &
Shipping

Product Images



Description

Whirl-Pak® bags with labeling field facilitate sample identification!

The white writing field can be labeled with an ordinary ballpoint pen or pencil.

Ideal for sampling liquids such as milk, water and food.

They can be frozen at any temperature, even in liquid nitrogen – careful handling is required after freezing.

Additional Information

No.	1-7163
Manufacturer (Brand)	Whirl Pak
EAN	4058072155797
Color	transparent (clear)
Material	Low-density polyethylene (LDPE)
Volume	710 ml
sterile	Yes
Bag features	with closure/lid with labeling field
Wide	150 mm
Length	230 mm
Type bag	Sample Bag Storage Bag
Type closure/ lid	Swing stopper
TBST MAX	82 °C
TBST MIN	-210 °C
Thickness	1.02 mm

