



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Whirl Pak® plastic pouch with title block, 125 x 305 mm (W x L), sterile

€78.40
plus VAT &
Shipping

Product Images



Description

Whirl-Pak® bags with labeling field facilitate sample identification!

The white writing field can be labeled with an ordinary ballpoint pen or pencil.

Ideal for sampling liquids such as milk, water and food.

They can be frozen at any temperature, even in liquid nitrogen – careful handling is required after freezing.

Additional Information

No.	1-7164
Manufacturer (Brand)	Whirl Pak
EAN	4058072154714
Color	transparent (clear)
Material	Low-density polyethylene (LDPE)
Volume	1260 ml
sterile	Yes
Bag features	with closure/lid with labeling field
Wide	125 mm
Length	305 mm
Type bag	Sample Bag Storage Bag
Type closure/ lid	Swing stopper
TBST MAX	82 °C
TBST MIN	-210 °C
Thickness	0.76 mm

