




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® storage box, 12 x 12 x 5.5 cm, 49 places

Special Price

€6.50 was

€18.60

**plus VAT &
Shipping**

Product Images



Description

Transparent storage box made of polystyrene with white plastic insert.

7 x 7 rows

Compartment size: 14 x 14 mm, suitable for vessels of 40–50 mm height, such as eppis, cryotubes 2 ml, microvials, etc.

Temperature resistant from -40°C to +70°C.

Stackable due to special shaping of the base and lid.

External dimensions: 12 x 12 x 5.5 cm high.

Additional Information

No.	2-1530
Manufacturer (Brand)	neoLab
EAN	4058072057282
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear) blue
Material	Polystyrene (PS)
Boxes properties	stackable/assemblable with removable grid insert with closure/lid with grid insert for different vessel types/volumes
Wide	120 mm
Height	55 mm
Length	120 mm
Depth	120 mm
Type closure/ lid	Hinged latch/lid
suitable for	Vials Micro Vials/ Tubes
Tray size	14 x 14 mm
Squares (WxL)	7 x 7
for vessels with DM up to	14.0 mm
for vessels with volume up to	2 ml
for type (micro) vessel	Cryotubes Reaction vessels
for type vial/vial	for threaded vials. 13 mm
suitable for cooling in	Freezer/ freezer

