



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® Storage box (PP) f. Microvials, 50 places, for 12 mm Ø

€7.96
plus VAT &
Shipping

Product Images



Description

Storage boxes with 5 x 10 spaces for vials with max. 12 mm Ø.
With organizer insert and tightly closing lid.

Suitable for e.g.:

Autosampler vials 1.2 ml (Order No. 2-7200 and 2-7201)

Micro vials 1.5 ml (Order No. 2-1000)

Conic ampoule made of PP (Order No. 2-1009)

Made of crystal clear polystyrene or translucent polypropylene, autoclavable.

Dimensions: 140 x 70 x 40 mm (L x W x H)

Additional Information

No.	2-1535
Manufacturer (Brand)	neoLab
EAN	4058072057299
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Boxes properties	with closure/lid for different vessel types/volumes
Wide	70 mm
Height	40 mm
Length	140 mm
Depth	140 mm
Type closure/ lid	Hinged latch/lid
suitable for	Vials/ Vials
for vessels with volume up to	1.2 ml 1.5 ml
suitable for cooling in	Freezer/ freezer
Tray size	12 x 12 mm
Squares (WxL)	5 x 10
for vessels with DM up to	12.0 mm
for type vial/vial	Vial. 12 mm Ø

