



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® Storage Box (PS) f. Microvials, 50 places, for 12 mm Ø

**€5.72**  
plus VAT &  
Shipping

### Product Images



## Description

---

Storage boxes with 5 x 10 spaces for vials with max. 12 mm Ø.  
With organizer insert and tightly closing lid.

Suitable for e.g.:

Autosampler vials 1.2 ml (Order No. 2-7200 and 2-7201)

Micro vials 1.5 ml (Order No. 2-1000)

Conic ampoule made of PP (Order No. 2-1009)

Made of crystal clear polystyrene or translucent polypropylene, autoclavable.

Dimensions: 140 x 70 x 40 mm (L x W x H)

## Additional Information

No.	2-7209
Manufacturer (Brand)	neoLab
EAN	4058072009038
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polystyrene (PS)
Boxes properties	with closure/lid with grid insert for different vessel types/volumes
Wide	70 mm
Height	40 mm
Length	140 mm
Depth	140 mm
Type closure/ lid	Hinged latch/lid
suitable for	Vials/ Vials
for vessels with volume up to	1.2 ml 1.5 ml
suitable for cooling in	Freezer/ freezer
Tray size	12 x 12 mm
Squares (WxL)	5 x 10
for vessels with DM up to	12.0 mm
for type vial/vial	Vial. 12 mm Ø

