



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton storage box, for 60 vials up to 13 mm Ø, 50 mm high

€36.71 As low as
€27.42
plus VAT &
Shipping

Product Images



Description

Compact storage box made of polystyrene with compartments for 40 and 60 counting vials, with foam insert at the bottom and hinged lid with alphanumeric table.

Order No. 9-0260

for 60 vials (max. 13 mm Ø, 50 mm high)

Accessories: Index cards (Order No. 9-0262 and 9-0263)

Additional Information

| | |
|-------------------------------|--|
| No. | WD-0856 |
| Manufacturer (Brand) | Wheaton |
| old neoLab article no. | 9-0260 |
| EAN | 4058072427245 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Color | transparent (clear) |
| Material | Polystyrene (PS) |
| Boxes properties | stackable/assemblable with removable grid insert with closure/lid with grid insert for different vessel types/volumes alphanumerically coded |
| Wide | 95 mm |
| Height | 60 mm |
| Length | 186 mm |
| Depth | 186 mm |
| Type closure/ lid | Hinged latch/lid |
| suitable for | Vials/ Vials |
| for vessels with volume up to | 2 ml 3 ml 4 ml |
| suitable for cooling in | Ultra Deep Freezer/ Deep Freezer |
| Tray size | 13 x 13 mm |
| Squares (WxL) | 6 x 10 |
| for vessels with DM up to | 13.0 mm |
| for type vial/vial | Vial. 12.5 mm Ø |

