



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine LaboBox system: box w. lid, PC, 39 mm inner height

**€12.10** As low as  
**€11.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

**The stackable LaboBox® boxes** with slip lid for holding the LaboBox® racks are available **in two different heights**.

Simply select according to the loading height of the rack.

The higher box is also suitable for slides 76 x 45 mm.

The boxes are **autoclavable and cryogenic**, temperature resistant from -130°C to +130°C.

**Material:** PC, transparent

**Dimensions without lid:** 170 x 83 mm (L x W)

**Order No. 2-2316**

Box with lid, inner height 39 mm

Dimensions: 178 x 83 x 43 mm (L x W x H)

## Additional Information

No.	2-2316
Manufacturer (Brand)	neoLabLine
EAN	4058072060541
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polycarbonate (PC)
Boxes properties	with closure/lid alphanumerically coded
Wide	83 mm
Width/inside	83 mm
Height	43 mm
Height inside	39.0 mm
Length	178 mm
Depth	178 mm
Depth inside	178 mm
Type closure/ lid	Hinged latch/lid
suitable for	Racks
suitable for cooling in	Freezer/ freezer
Tray size	2 x 76 mm
Squares (WxL)	1 x 50
for type (micro) vessel	Reaction vessels Cryogenic tubes

