



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Mini-Box transparent, PE, 350 ml, 130 x
60 x 60 mm**

**€2.50
plus VAT &
Shipping**

Product Images



Description

Storageboxes with rounded corners.

- sturdy design, stackable
- tightly closing snap-on lid, with labeling field
- completely transparent

Material: polyethylene

Additional Information

No.	2-1684
Manufacturer (Brand)	neoLab
EAN	4058072072421
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	350 ml
sterile	No
Boxes properties	with closure/lid
Wide	60 mm
Form	rectangular
Height	60 mm
Length	130 mm
Depth	130 mm
Type closure/ lid	Push-on/plug-in closure
suitable for	variable applications
TBST MAX	80 °C
TBST MIN	-40 °C
Closure / lid properties	closed
suitable for cooling in	Freezer/ freezer
Color closure/ lid	blue transparent (milky)
Material closure/ lid	Low-density polyethylene (LDPE)

