




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Multifunction box large, 1.7 L, transparent

€8.47
plus VAT &
Shipping

Product Images



Description

The neoLab multifunctional boxes made of robust PP in a beautiful design are the little helpers in the laboratory. The boxes with attached hinged lid are suitable for storing and transporting various samples, small parts, tools and other utensils needed in everyday laboratory work.

- dishwasher safe

Additional Information

No.	2-1527
Manufacturer (Brand)	neoLab
EAN	4058072118686
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	1700 ml
sterile	No
Boxes properties	with closure/lid
Wide	145 mm
Form	rectangular
Height	65 mm
Length	245 mm
Depth	245 mm
Type closure/ lid	Hinged latch/lid
Closure / lid properties	closed
suitable for cooling in	Refrigerator freezer/ freezer
Color closure/ lid	transparent (milky)
Material closure/ lid	Polypropylene (PP)

