




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Multifunction box small, 0.9 L, green

€6.05
plus VAT &
Shipping

Product Images

Description

The neoLab multifunctional boxes made of robust PP in a beautiful design are the little helpers in the laboratory. The boxes with attached hinged lid are suitable for storing and transporting various samples, small parts, tools and other utensils needed in everyday laboratory work.

- dishwasher safe

Additional Information

No.	2-1537
Manufacturer (Brand)	neoLab
EAN	4058072118662
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (light green)
Material	Polypropylene (PP)
Volume	900 ml
sterile	No
Boxes properties	with closure/lid
Wide	129 mm
Form	rectangular
Height	60 mm
Length	177 mm
Depth	177 mm
Type closure/ lid	Hinged latch/lid
Color closure/ lid	green
Material closure/ lid	Polypropylene (PP)
Closure / lid properties	closed
suitable for cooling in	Refrigerator freezer/ freezer

