



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® freezer can made of PP, 1.7 l

€6.08 As low as
€5.60
plus VAT &
Shipping

Product Images



Description

Square and rectangular all-purpose jars for storage and freezing with tight-fitting elastic lid.
The jars are transparent with blue lid.
Stackable.

Material: polypropylene

Additional Information

No.	2-1587
Manufacturer (Brand)	neoLab
EAN	4058072057534
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	1700 ml
sterile	No
Boxes properties	with closure/lid
Wide	120 mm
Form	rectangular
Height	80 mm
Length	250 mm
Depth	250 mm
Type closure/ lid	Push-on/plug-in closure
suitable for	variable applications
Color closure/ lid	blue
Material closure/ lid	Polypropylene (PP)
TBST MAX	110 °C
TBST MIN	-40 °C
Closure / lid properties	closed
suitable for cooling in	Freezer/ freezer

