



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cryo Transport Box Styrofoam large, 375 x 205 x 145 mm, 10 pcs./pack


Special Price

€13.30 was

€19.95

**plus VAT &
Shipping**

Product Images



Description

Lightweight, very stable cryo and transport box for tubes and vials.
Styrofoam 100% CFC-free .

Dimensions (WxDxH):
355x180x70 mm (inside)
375x205x90 mm (outside)

Packaging unit: 10 pieces/pack

Additional Information

No.	4-5021
Manufacturer (Brand)	neoLab
EAN	4058072023492
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Polystyrene (EPS. Styropor®)
Boxes properties	with removable grid insert with closure/lid for light-sensitive specimens
Wide	205 mm
Width/inside	180 mm
Height	90 mm
Height inside	70.0 mm
Length	375 mm
Depth	375 mm
Depth inside	355 mm
Type closure/ lid	Slip closure
for vessels with volume up to	50 ml
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
for vessels with DM up to	30.0 mm

