




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine storage box for 50 ml tubes, 3 x 3 places, orange

€14.80
plus VAT &
Shipping

Product Images



Description

Storage box for 15/50 ml cell culture tubes (Falcon) in bright colors.

Two inserts ensure a firm hold of the tubes in the upper as well as in the lower area.
The boxes are temperature resistant down to -80°C and have a transparent, milky lid.
5 x 5 places for 15 ml cell culture tubes with max. 17 mm diameter.
3 x 3 places for 50 ml cell culture tubes with max. 30 mm diameter.

Additional Information

No.	2-1909
Manufacturer (Brand)	neoLabLine
EAN	4058072058784
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	orange
Material	Polypropylene (PP)
Boxes properties	stackable/ collapsible with removable louvre insert with closure/lid with louvre insert
Wide	137 mm
Height	128 mm
Length	137 mm
Depth	137 mm
Type closure/ lid	Hinged latch/lid
suitable for	Centrifuge tube
Capacity	9
Squares (WxL)	3 x 3
for vessels with DM up to	30.0 mm
for vessels with volume up to	50 ml
for type (micro) vessel	Centrifuge tube
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

