



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine storage box for 15 ml tubes, 5 x 5 places, transparent

€12.30
plus VAT &
Shipping

Product Images



Description

Storage box for 15/50 ml cell culture tubes (Falcon) in bright colors.

Two inserts ensure a firm hold of the tubes in the upper as well as in the lower area.

The boxes are temperature resistant down to -80°C and have a transparent, milky lid.

5 x 5 places for 15 ml cell culture tubes with max. 17 mm diameter.

3 x 3 places for 50 ml cell culture tubes with max. 30 mm diameter.

Additional Information

No.	2-1917
Manufacturer (Brand)	neoLabLine
EAN	4058072058869
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Boxes properties	stackable/ collapsible with removable louvre insert with closure/lid with louvre insert
Wide	137 mm
Height	128 mm
Length	137 mm
Depth	137 mm
Type closure/ lid	Hinged latch/lid
suitable for	Centrifuge tube
Squares (WxL)	5 x 5
for vessels with DM up to	17.0 mm
for vessels with volume up to	15 ml
for type (micro) vessel	Centrifuge tube
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

